#### MCE USA

Tel: 609-228-6898

Email: tech@MedChemExpress.com



# JC-1 Mitochondrial Membrane Potential Assay Kit

## 1 Contents

Components	JC-1(200 μM in DMSO)	Phosphate-buffered saline (10×)	CCCP (50 mM in DMSO)
HY-K0601-100T	230 μL × 5	25 mL	125 μL

#### 2 General Information

MCE JC-1 Mitochondrial Membrane Potential Assay Kit uses JC-1, a lipophilic cationic dye, to detect the mitochondrial membrane potential in variety of cell types, as well as intact tissues and isolated mitochondria.

 $\Delta\Psi$ m, mitochondrial membrane potential, is an important parameter of mitochondrial function and has been used as an indicator of cell health. Variation of  $\Delta\Psi$  m would be studied using JC-1. In healthy cells with high  $\Delta\Psi$ m, JC-1 forms complexes known as J-aggregates. While in cells with low  $\Delta\Psi$ m, JC-1 remains in the monomeric form. When excited at 490 nm, JC-1 monomers emit a green fluorescence with a maximum at ~520 nm. Aggregates of JC-1 emit an orange-red fluorescence with a maximum at ~590 nm.

Therefore, the fluorescence intensity of the orange-red emission and the ratio of orange-red fluorescence to green fluorescence can be used to measure mitochondrial membrane potential and serve as an indicator of overall cell health.

MCE JC-1 Mitochondrial Membrane Potential Assay Kit consists of JC-1, Phosphate-buffered saline (10×) and CCCP (50 mM in DMSO). CCCP controls should be used to confirm that the JC-1 response is sensitive to changes in membrane potential.

MCE JC-1 Mitochondrial Membrane Potential Assay Kit can be analyzed by fluorescence microscopy, flow cytometer, or a fluorescence plate reader with appropriate filter sets.

#### 3 General Protocol

- 1. Preparation of Phosphate-buffered saline (PBS) (1×)
- a. Warm the PBS (10×) until any salt crystals are completely dissolved.
- b. Dilute PBS (10×) 1:10 with dH<sub>2</sub>O (e.g. 1 mL PBS + 9 mL dH<sub>2</sub>O).
- 2. Labeling of Cells
- a. Culture cells in 6-, 12-, 24-, or 96-well plates at a density of 5-10× 10<sup>s</sup> cells/mL. Incubate the cells according to your normal protocol.
- b. For the control tube, allow the vial of CCCP has come to room temperature, add 1 μL of CCCP (50 mM). Incubate cells at 37 °C for 5 minutes.
- c. Add 10 μL JC-1 (200 μM) per well to make the final concentration at 2 μM. Incubate cells at 37 °C, 5% CO<sub>2</sub>, for 15-20 minutes.

Note: If additional labeling followed, for example with an annexin V, begin with step 3.a.

- d. After incubation, centrifuge cells for 3-4 minutes at 400× g at 4°C, carefully aspirate the supernant.
- e. Wash cells twice with PBS (1 $\times$ ): add 2 mL PBS (1 $\times$ ) to suspend cells and vortex to mix thoroughly. Centrifuge cells for 3-4 minutes at 400 $\times$  g at 4 $^{\circ}$ C, carefully aspirate the supernant.
- f. Add 500 µL PBS (1×) to suspend cells. Analyze sample on a flow cytometer, fluorescence microscopy, or fluorescence microplate reader.
- 3. Additional Labeling with Annexin V
- a. After step 2.c, wash cells twice with PBS (1×): add 2 mL PBS (1×) to suspend cells and vortex to mix thoroughly. Centrifuge cells for 3-4 minutes at 400× g at
- $4\,^\circ\!\!\mathrm{C}$  , carefully aspirate the supernant.
- b. Add 100 µL Annexin V binding buffer (10 mM HEPES, 140 mM NaCl, 2.5 mM CaCl<sub>2</sub>, pH 7.4) to suspend the JC-1-stained cells.
- c. Add appropriate volume of Annexin V and incubate cells at 37°C, 5% CO<sub>2</sub>, for 15 minutes.

d. Add 400 μL Annexin V binding buffer (10 mM HEPES, 140 mM NaCl, 2.5 mM CaCl<sub>2</sub>, pH 7.4) to suspend cells. Analyze sample on a flow cytometer, fluorescence microscopy, or fluorescence microplate reader.

#### 4 Storage

Stored at -20  $^{\circ}$  protecting from light, and is stable for up to 12 months.

For immediate use, components may be stored at 2-8°C.

### 5 Precautions

- 1. A 'Test' simply refers to a single assay well. 100 tests based on labeling volumes of 1.0 mL.
- 2. JC-1 is light sensitive, we recommend that you make small aliquots and store them at -20 C with light protected.
- 3. Phosphate-buffered saline (10×) could be stored at 4°C.
- 4. This product is for R&D use only, not for drug, house hold, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Tel: 609-228-6898