

# PolyFast Transfection Reagent

## 1 Contents

Components	HY-K1014-750 $\mu$ L	HY-K1014-1.5 mL	HY-K1014-7.5 mL
PolyFast Transfection Reagent	750 $\mu$ L	1.5 mL	5 x 1.5 mL

## 2 Introduction

MCE PolyFast Transfection Reagent consists of cationic polymers and can introduce nucleic acids (DNA or RNA) into eukaryotic cells. In polymer-based transfection, exogenous DNA forms complexes with cationic polymers that enter host cells by endocytosis. PolyFast Transfection Reagent has high-efficiency, low-toxicity transfection of many cell types, including some hard-to-transfect cells.

## 3 General Protocol

### 1. Prepare cells

1.1 Inoculate cells in advance until the density reaches 70-90% for cell transfection.

Note: The viability and general health of cells prior to transfection significantly affect transfection result. Cells should be at least 90% viable prior to transfection and have had sufficient time to recover from passaging.

1.2 Remove the medium. Wash twice with PBS, and then add 900  $\mu$ L serum-free medium to each well of 6-well plate.

### 2. Prepare PolyFast Transfection Reagent/DNA complex

2.1 According to table 1, dilute 3  $\mu$ L PolyFast Transfection Reagent with 50  $\mu$ L serum-free medium for each well of 6-well plate and mix gently. Incubate at room temperature for 5 minutes.

2.2 According to table 1, dilute 1  $\mu$ g DNA with 50  $\mu$ L serum-free medium for each well of 6-well plate and mix gently. Incubate at room temperature for 5 minutes.

2.3 Mix the diluted PolyFast Transfection Reagent and DNA gently. Incubate at room temperature for 15 minutes.

### 3. Add PolyFast Transfection Reagent/DNA complex to cells

Add the PolyFast Transfection Reagent/DNA complex to cells in 6-well plate and mix well. Incubate at 37°C for 24-48 hours. The medium can be replaced with fresh serum-containing medium after 6 hours if necessary.

Table 1. PolyFast Transfection Reagent: DNA Ratio.

Plate Size		96-well	24-well	12-well	6-well	60 mm	100 mm
Growth Area (per well)		0.3 cm <sup>2</sup>	2 cm <sup>2</sup>	4 cm <sup>2</sup>	9.5 cm <sup>2</sup>	20 cm <sup>2</sup>	60 cm <sup>2</sup>
Serum-free medium		80 $\mu$ L	450 $\mu$ L	630 $\mu$ L	900 $\mu$ L	2.7 mL	5.4 mL
Dilution of Transfection Reagent	Transfection Reagent	0.75 $\mu$ L	1.5 $\mu$ L	2.25 $\mu$ L	3 $\mu$ L	7.5 $\mu$ L	15 $\mu$ L
	Serum-free medium	10 $\mu$ L	25 $\mu$ L	35 $\mu$ L	50 $\mu$ L	150 $\mu$ L	300 $\mu$ L
Dilution of DNA or RNA	DNA	0.25 $\mu$ g	0.5 $\mu$ g	0.75 $\mu$ g	1 $\mu$ g	2.5 $\mu$ g	5 $\mu$ g
	Serum-free medium	10 $\mu$ L	25 $\mu$ L	35 $\mu$ L	50 $\mu$ L	150 $\mu$ L	300 $\mu$ L
Total Volume (per well)		100 $\mu$ L	500 $\mu$ L	700 $\mu$ L	1 mL	3 mL	6 mL

Note: 1). Generally, the ratio of DNA ( $\mu\text{g}$ ) to PolyFast Transfection Reagent ( $\mu\text{L}$ ) is 1:3, and the transfection effect can be optimized in the range of 1:1 to 1:5 if necessary.

2) RNA transfection just follows the protocol as described for DNA.

#### 4. Detection

Measure transfection efficiency using an appropriate assay for the reporter gene. For transient transfection, cells are typically assayed 24-48 hours after transfection.

#### 4 Storage

Store at  $-20^{\circ}\text{C}$  for 2 years. Avoid repetitive freeze-thaw cycles.

#### 5 Precautions

1. The viability and general health of cells prior to transfection significantly affect transfection result. Cells should be at least 90% viable prior to transfection and have had sufficient time to recover from passaging.
2. Take high-purity, sterile, contaminant-free DNA or RNA for transfection.
3. Dilute the PolyFast Transfection Reagent and DNA in serum-free medium.
4. This product is for R&D use only, not for drug, household, or other uses.
5. For your safety and health, please wear a lab coat and disposable gloves to operate.