


LY2109761

Cat. No.: HY-12075
Batch No.: 32412
CAS No.: 700874-71-1
Molecular Formula: C₂₆H₂₇N₅O₂

ANALYTICAL DATA

Theoretical: C(70.73%), H(6.16%), N(15.86%)
Result: C(70.71%), H(6.20%), N(15.75%)

Chemical Structure:



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA