

Product Data Sheet

Biotin-probe 1

Cat. No.: HY-135641 Molecular Formula: $C_{32}H_{58}N_8O_8S$ Molecular Weight: 714.92

Target: Fluorescent Dye

Pathway: Others

Storage: Please store the product under the recommended conditions in the Certificate of

Analysis.



BIOLOGICAL ACTIVITY

Description

Biotin-probe 1 is a non-radiolabeled probe. Biotin-labeled probes can be applied to in situ hybridization^[1].

REFERENCES

[1]. Bloch B.Biotinylated probes for in situ hybridization histochemistry: use for mRNA detection. J Histochem Cytochem. 1993 Dec;41(12):1751-4.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA

Screening Libraries

Inhibitors

יי