# **Product** Data Sheet

## Cloxacillin

**Cat. No.:** HY-B0466A **CAS No.:** 61-72-3

Molecular Weight: 435.88

Target: Antibiotic; Bacterial; Beta-lactamase

Pathway: Anti-infection

**Storage:** Please store the product under the recommended conditions in the Certificate of

Analysis.

## **BIOLOGICAL ACTIVITY**

Description	Cloxacillin is an orally active antibacterial agent and $\beta$ -lactamase inhibitor with an IC <sub>50</sub> of 0.04 $\mu$ M. Cloxacillin can suppress the <b>S.</b> <i>aureus</i> -induced inflammatory response by inhibiting the activation of MAPKs, NF- $\kappa$ B and NLRP3-related proteins <sup>[1][2]</sup> [3].				
IC <sub>50</sub> & Target	Tetracycline				
In Vitro	Cloxacillin (0-2048 $\mu$ g/mL; 20-24 h) shows good antibacterial activity for S. aureus 8325-4 and DU1090 with MIC values both of 0.125 $\mu$ g/mL <sup>[1]</sup> . Cloxacillin (0.015625 $\mu$ g/mL; 6 h) inhibits the hemolytic activity of HI $\alpha$ in vitro, and this inhibition is not only more pronounced when combined with TZ and TZ, but also suppresses the inflammatory response by inhibiting the activation of MAPKs, NF- $\kappa$ B and NLRP3-related proteins <sup>[1]</sup> . MCE has not independently confirmed the accuracy of these methods. They are for reference only. Cell Viability Assay <sup>[1]</sup>				
	Cell Line:	S. aureus 8325-4, S. aureus DU1090 (an Hlα-deleted strain)			
	Concentration:	0-2048 μg/mL			
	Incubation Time:	20-24 h			
	Result:	Inhibited S. aureus 8325-4 and DU1090 with MIC values both of 0.125 $\mu g/mL$ .			
	Western Blot Analysis <sup>[1]</sup>				
	Cell Line:	S. aureus 8325-4			
	Concentration:	0.015625 μg/mL (combines with Thioridazine (TZ, 0.25 μg/mL) and Tetracycline (TC, 0.03125 μg/mL)).			
	Incubation Time:	6 h			
	Result:	Inhibited the expression of Hl $\!\alpha$ and the inhibition was more pronounced when combined with TZ and TC.			
	Western Blot Analysis <sup>[1]</sup>				

Cell Line:	RAW264.7 cells (exposes to S. aureus 8325-4/DU1090 or pure Hlα)
Concentration:	$0.015625\mu g/mL$ (combines with TZ (0.25 $\mu g/mL)$ and TC (0.03125 $\mu g/mL)).$
Incubation Time:	6 h
Result:	Inhibited the activation of MAPKs, NF-κB and NLRP3-related proteins thereby inhibiting the inflammatory response when combined with TC and TZ.

#### In Vivo

Cloxacillin (1.6125 mg/kg; s.c.; 12-h intervals for 72 h) protects mice from S. aureus peritonitis in vivo when combines with Thioridazine and Tetracycline  $^{[1]}$ .

Cloxacillin (7.5 mg/per; i.p.; twice daily from day 3 for 3 days) develops less severe synovitis and reduces bone erosions when combines with anti-IL-15 antibodies  $^{[3]}$ .

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

Female BALB/c mice (6-week-old; peritonitis model) <sup>[1]</sup> .			
1.6125 mg/kg (combines with TC (3.125 mg/kg) and TZ (25 mg/kg))			
Subcutaneous injection; 12-h intervals for 72 h.			
Reduced the degree of inflammatory cell infiltration in the mouse lung tissue and alveolar structures tended to be normal.			
Significantly reduced the pathological changes in spleen and liver tissue, as well as decreased the CFU counts of S. aureus in the peritoneal cavity.			

Animal Model:	Female wildtype C57BL/6 mice (8-week-old; systemic S. aureus-induced arthritis model)
Dosage:	7.5 mg/per (combines with 25 μg/per anti-IL-15 antibodies)
Administration:	Intraperitoneal injection; twice daily from day 3 (after bacterial inoculation) and stopped at day 6.
Result:	Showed activities of reducing severe synovitis and bone erosions when combined with anti-IL-15 antibodies.

## **CUSTOMER VALIDATION**

- Chemosphere. 2019 Jun;225:378-387.
- Int J Biol Macromol. 25 December 2021.
- Biomed Res Int. 2018 Jul 2;2018:3579832.

See more customer validations on www.MedChemExpress.com

### **REFERENCES**

[1]. Zhou H, et al. The combination of cloxacillin, thioridazine and tetracycline protects mice against Staphylococcus aureus peritonitis by inhibiting  $\alpha$ -Hemolysin-induced MAPK/NF- $\kappa$ B/NLRP3 activation. Int J Biol Macromol. 2022 Feb 15;198:1-10.

[2]. Bergmann B, et al. Antibiotics with Interleukin-15 Inhibition Reduce Joint Inflammation and Bone Erosions but Not Cartilage Destruction in Staphylococcus aureus-

Page 2 of 3 www.MedChemExpress.com

Induced Arthritis. Infect Immun. 2018 Apr 23;86(5):e00960-17.						
[3]. Lupiola-Gómez PA, et al. Gro	oup 1 beta-lactamases of Aeron	nonas caviae and their resistance	to beta-lactam antibiotics. Can J Micr	robiol. 2003 Mar;49(3):207-15.		
			al applications. For research use			
	Tel: 609-228-6898 Address: 1 De	Fax: 609-228-5909 er Park Dr, Suite Q, Monmouth	E-mail: tech@MedChemExpress.n Junction, NJ 08852, USA	.com		

Page 3 of 3 www.MedChemExpress.com