


Carbostyryl 165

Cat. No.:	HY-D1572
CAS No.:	26078-23-9
Molecular Formula:	C ₁₂ H ₁₄ N ₂ O
Molecular Weight:	202.25
Target:	Fluorescent Dye
Pathway:	Others
Storage:	Please store the product under the recommended conditions in the Certificate of Analysis.



BIOLOGICAL ACTIVITY

Description

Carbostyryl 165 is an electronically neutral fluorescent dye (Ex: 351/364 nm, Em: >400 nm)^[1].

REFERENCES

[1]. Windsor S A, et al. Electro-fluorescence polarization studies of the interaction of fluorescent dyes with clay minerals in suspensions. *Colloids and Surfaces A: Physicochemical and Engineering Aspects*, 1999, 148(1-2): 61-73.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA