


## Sulfo-Cy3.5 maleimide

Cat. No.:	HY-D2239
Molecular Formula:	C <sub>44</sub> H <sub>43</sub> K <sub>3</sub> N <sub>4</sub> O <sub>15</sub> S <sub>4</sub>
Molecular Weight:	1113.38
Target:	Fluorescent Dye
Pathway:	Others
Storage:	Please store the product under the recommended conditions in the Certificate of Analysis.



### BIOLOGICAL ACTIVITY

#### Description

Sulfo-Cy3.5 maleimide is sulfonated modified Cy3.5 with good water solubility.  $\lambda_{ex}$ =591 nm,  $\lambda_{em}$ =604 nm<sup>[1]</sup>.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA