

Cy5- β -Amyloid (42-1), human

Cat. No.:	HY-P1362F	
Target:	Amyloid- β	
Pathway:	Neuronal Signaling	
Storage:	-20°C, sealed storage, away from moisture and light * In solvent : -80°C, 6 months; -20°C, 1 month (sealed storage, away from moisture and light)	Cy5- β -Amyloid (42-1), human

SOLVENT & SOLUBILITY

In Vitro	H ₂ O : 20 mg/mL (Need ultrasonic)
----------	---

BIOLOGICAL ACTIVITY

Description	Cy5- β -Amyloid (42-1), human is a Cy5 fluorescently-labelled β -Amyloid (42-1, human) peptide (λ_{ex} = 633 nm and λ_{em} = 670 nm) ^[1] .
-------------	--

REFERENCES

[1]. Jungbauer LM, et al. Preparation of fluorescently-labeled amyloid-beta peptide assemblies: the effect of fluorophore conjugation on structure and function. J Mol Recognit. 2009 Sep-Oct;22(5):403-13.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA