

## CDKN2A/p14ARF Antibody (YA1563)

Cat. No.:	HY-P81818
Host:	Rabbit
Reactivity:	Human
Conjugation:	Non-conjugated
SwissProt ID:	Q8N726
Research Field:	Cell Biology
Molecular Weight:	Predicted band size: 14 kDa; Observed band size: 14 kDa

### PROPERTIES

Formulation	Supplied in rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.	
Purity	Affinity Chromatography	
Storage & Stability	Stored at -20°C for 1 year. Avoid repeated freeze / thaw cycles.	
Appearance	Liquid	
Application & Dilution Ratio	Application	Dilution Ratio
	WB	1/500-1/1000
	IF	1/50-1/200
	IP	1/50
	FC	1/50-1/100
Shipping	Shipping with blue ice.	

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA