

IKB beta Antibody (YA3088)

Cat. No.:	HY-P83343
Host:	Rabbit
Reactivity:	Human, Mouse
Conjugation:	Non-conjugated
SwissProt ID:	Q15653
Research Field:	Signal Transduction
Molecular Weight:	Predicted band size: 38 kDa; Observed band size: 48 kDa

PROPERTIES

Formulation	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA	
Purity	Affinity Purified	
Storage & Stability	Stored at -20°C for 1 year. Avoid repeated freeze / thaw cycles.	
Appearance	Liquid	
Application & Dilution Ratio	Application	Dilution Ratio
	WB	1/500-1/1000
	IHC	1/50-1/100
	IP	1/20
Shipping	Shipping with blue ice.	

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA